


Operator's Manual

Spectroil M/N-W & Spectroil M/C-W Oil Analysis Spectrometers

(For Serial Numbers Starting with 6001, Modifications 6 & 7)



SPECTROINC.

QinetiQ North America

160 Ayer Road Littleton, MA 01460 U.S.A. • Tel. 978.486.0123 • Fax 978.486.0030
email: service@spectroinc.com • web page: www.spectroinc.com

Table of Contents

1.0	GENERAL OPERATING REQUIREMENTS	1
1.1	Power Application and Systematic Power Removal.....	1
1.2	Rod Electrode Sharpening	3
1.3	Installing the Disc Electrode.....	3
1.4	Installing the Rod Electrode and Setting the Gap	4
1.5	Installing and Positioning the Sample Holder	5
1.5.1	Disposable Sample Holder	5
1.5.2	Reusable Sample Holder.....	6
1.5.3	Sample Holder Cover	6
1.6	Cleaning the Sample Stand	7
1.6.1	Cleaning After Each Burn Cycle	7
1.6.2	Cleaning After Each Operating Shift	7
1.6.3	Cleaning the Quartz Window	7
1.6.4	Cleaning Solutions	8
1.7	Paper Tissue for Operating and Paper Towels for Cleaning	8
1.8	Waste Oil Disposal Container	8
2.0	DAILY OPERATION	8
2.1	Daily Routine Prior to Use	8
2.2	Warm-Up Procedure	10
2.3	Routine Sample Analysis	10
2.4	Daily Standardization Check	11
2.5	Complete Standardization	12
2.6	Daily Routine Prior to Securing	16
2.7	Optical Profiling	17
2.8	Sample Identification (I.D.)	17
2.8.1	Setup of the Sample I.D.	18
2.8.2	Using Sample ID.'s	18
2.9	Sulfur Analysis.....	19
2.9.1	Set-up for Sulfur Analysis.....	19
2.9.2	Profiling the sulfur Optic	21
2.9.3	Standardization of the Sulfur Optic.....	21
2.9.3	Verification	21
2.9.4	Clean Up and Final Checks	22
2.10	Coolant Analysis	22
2.10.1	Coolant Program Standardization	22
2.10.2	Verification	24
2.10.3	Routine Coolant Sample Analysis.....	24
2.10.4	Grease Analysis.....	24
3.0	OPERATOR MAINTENANCE	25
3.1	Daily Operator Maintenance.....	25
3.2	Scheduled Periodic Maintenance	25
3.3	Procedure to Replace Electrode Sharpener Cutter Blade, M90102.....	30

4.0	PERFORMING CALIBRATION CURVE VERIFICATION	30
4.1	Repeatability Testing	30
4.1.1	Repeatability Specifications	31
4.1.2	Repeatability Test	31
4.1.3	Factors Affecting Repeatability	32
4.2	Accuracy Testing	32
4.3	Disc Electrode Offset Procedure	33
5.0	BACKUP OF COMPLETE OILM SOFTWARE	35

List of Effective Pages

Summary of Spectroil Modifications


Modification	Description
Mod 0	Original CID version of Spectroil M
Mod 1	Upgrade with SFTM port and frequency adjustment potentiometer
Mod 2	Addition of solid state excitation ignition module and SFTM port
Mod 3	Upgrade to combined solid state source
OilMWindows	Upgrade to Windows hardware and software
Mod 4	Upgrade to panel PC hardware and Windows XP
Mod 5	CE version of the Spectroil M
Mod 6	CCD Optic and New Software v. 5, starting w/serial number 6001
Mod 7	C.E. version of Mod. 6

Summary of OilMWindows Modification 6 Hardware and Software Manual Versions

Change	Version	Date	Description
First Issue	3.0	7/1/07	Complete update of Mod. 5 version 2.5 manual to include CCD optic & updated software.
Change 1	3.1	10/24/07	Update for CE version, Mod. 7, changes. Addition of optic removal procedure. Grammatical corrections
Change 2	3.2	10/31/07	More updates to Chapter 2 for CE certification.
Change 3	3.3	2/18/08	Created separate Operator's Manual, updated accuracy and repeatability tables with additional elements, included sample ID and operator's maintenance.
Change 4	3.4	4/14/08	Added coolant analysis procedure; minor corrections
Change 5	3.4.1	4/15/08	Corrections to coolant analysis and sulfur procedure.
Change 6	3.4.2	7/21/09	Added grease analysis and minor corrections.
Change 7	3.4.3	12/23/09	Corrected Fig. 2-1 & replaced sample cap with sample holder.

Total Number of pages in this manual is 44 consisting of the following:

Section	Page Numbers
Cover	2
Table of Contents	i - ii
List of Effective Pages	iii - vi
Chapters 1 - 5	1 - 36



WARNING!!!

High Voltages are Present During the Operation of the Spectroil M!

Observe all Safety Precautions!

Turn OFF the Main Power Switch and unplug the SPECTROIL M before any work is performed.

Definitions

The following definitions apply to specific instructions throughout this manual.



WARNING!!!

An operating procedure or practice that may cause injury if not carefully observed or followed.

CAUTION!!!

An operating procedure or practice that may cause damage to the Spectroil M if not carefully observed or followed

NOTE!!!

An operating procedure or practice that is essential to emphasize

Software CAUTION!!!

The Spectroil M computer is capable of running multiple software applications and/or operating systems. However, as designed, the computer processor is dedicated to the operation and control functions of the Spectroil M. Do not attempt to add any software or alter the original factory installed software without checking first with the Spectro Inc. Service Department.

Note on Oil Standards

The Spectroil M series of spectrometers can be calibrated for military or commercial applications. As a rule, the Spectroil M/N-W is calibrated and standardized with D-19, D-12 and D-3 series of standards, and the Spectroil M/C-W with V-21 or S-21 series of calibration standards.

Although this manual frequently refers to the military “D” series of standards, the operator procedures are identical for all types of Spectroil M spectrometers. Commercial customers should substitute their equivalent “V” or “S” series of standards throughout this manual.

WARRANTY

The warranty period of the Spectroil M family of spectrometers is twelve (12) months from date of installation or fifteen (15) months from date of shipment, whichever occurs first. Spectro warrants the Spectroil under conditions of operation against defects of materials and workmanship. All defective material will be replaced providing damage was not caused by improper use. Warranty applies to parts and labor only.

Spectroil M/C-W

&

Spectroil M/N-W

Operator's Manual

1.0 GENERAL OPERATING REQUIREMENTS

This manual provides the Spectroil M operator with routine instructions on how to set-up the spectrometer, analyze samples, and perform routine maintenance. The instructions apply to the Spectroil M/C-W and the Spectroil M/N-W Oil Analysis Spectrometers with Mod. 6 (CCD optic) and Mod. 7 (C.E. hardware).

1.1 Power Application and Systematic Power Removal

The Spectroil M consists of the excitation, optics and readout subassemblies. Each subassembly requires specific voltages to perform as a system. These voltages are generated within each subassembly and originate from the main power distribution assembly. Main power for these subassemblies is fused with a 10 ampere circuit breaker CB1 which is mounted on the Power Connection Plate located on the right side near the back of the instrument. Power should not be applied to the instrument unless all specifications for input power requirements have been met.

To apply power to the instrument, place circuit breaker CB1 in the upright ON position, Figure 1-1.

Once power is applied to the instrument, two events occur. First, an internal controller will




Figure 1-1, Right Side View Showing Circuit Breaker CB1

initialize and automatically load the application. Simultaneously, the readout system will load Windows XP® and start up the instrument's application program called OILMWindows®. The instrument will boot directly to the Analysis Screen, Figure 1-2.

If the system fails to establish communication, a screen similar to Figure 1-3 appears to select the Correct Configuration File Path. This will happen if the system cannot find a file in case it has been moved, updated or is corrupted. Refer to the Maintenance Manual Section 2.4.4.6 for assistance to diagnose and correct this condition.

Next, move the MODE switch to the OPERATE position and power will be applied to the excitation source and electrode sharpener. A noticeable increase in fan activity will be observed; this is normal.

At this point, the instrument is ready to begin operation. Some time will be needed before the instrument stabilizes after power is applied. It is recommended that the main power remain on when the instrument is not in use to maintain




Figure 1-2, Analysis Screen




Figure 1-3, Error Message if Communications Cannot be Established

maximum instrument stability. When not in use, the MODE switch should be placed in the STANDBY position.

To completely turn off the Spectroil M, first shut down the OILMWindows software and then remove power from the instrument.

CAUTION: Utilizing the Windows XP® operating system dictates a specific series of steps to be performed in the process of shutting down the Spectroil M. If power is accidentally removed from the Spectroil M spectrometer or the circuit breaker CB1 is shut off while the readout system is running Windows XP® or the OILMWindows® application, an orderly shutdown would not be performed and as a result, the Windows XP® operating system must perform hardware and software diagnostics when power is reapplied.

To prevent this from happening, always follow the next steps to shut down the Spectroil M.

1. Shut down the OILMWindows® application by left clicking the close box (box with an "X") in the upper right corner of the OILMWindows® header, or choose File/Exit from the pull down menu options. This will return the software to the Windows XP® desktop.
2. Left click the START menu and select Shut Down. A dialog will appear asking what do you what the computer to do. The options are, Standby, Restart, or Shutdown.
3. Highlight the Shutdown radio button and select OK.
4. The Windows XP® logo will appear with the instruction that it is shutting down.
5. After the logo disappears a message "It's now safe to turn off your computer" will appear.
6. It is now safe to place the main circuit breaker CB1 in the down position to remove all

power from the instrument.

If, by accident, power was removed from the instrument, a series of diagnostics will automatically be performed when power is reapplied. It is extremely important to allow these diagnostics to complete in their entirety before loading any application software. If problems are experienced during the process of running these diagnostics, contact the technical service department of Spectro Incorporated for instructions on how to recover and proceed with normal operation.

To restore power to the instrument, follow the instructions given above.

1.2 Rod Electrode Sharpening

The rod electrode, along with the disc electrode, form the analytical gap through which the oil or fuel sample are passed for analysis. An alternating current discharge will occur between the disc and rod electrode and vaporize the sample and the metallic components in it. This is the basis of the arc emission technique.

The preparation of the tip of the rod electrode plays a significant role in obtaining repeatable analytical data. The rod electrode must be cleaned prior to inserting it into the electrode sharpener. This is accomplished by taking a clean paper towel and removing the components of the burn residue from the previous analysis. Remove all residue from the tip and sides of the electrode by rotating the rod in the paper towel while applying pressure with the fingers of the opposite hand.

NOTE: The paper towel should be laboratory grade and free of silicon.

With the spectrometer on, turn the MODE switch to OPERATE. This applies power to J2, the electrode sharpener power connector. Momentarily press the power switch located on the base of the electrode sharpener. The sound of the electric motor should be heard and a slight vibration should

be felt through the motor. The motor will continue to run on a self-timed cycle for approximately 3 to 5 minutes. To sharpen the rod electrode, insert the rod into the rotating electrode guide hole until it comes in contact with the cutter blade. Apply inward pressure until approximately 1/8 to 3/16 inch (3 to 5 mm) is cut from the end of the rod. Slightly decrease the inward pressure on the rod electrode, but still maintain its contact with the cutter blade. This will polish the rod electrode tip.

Remove the rod electrode and visually inspect the tip. It should have a clean cut with no apparent chipping around the circumference of the rod. The surface should be very smooth and have a polished mirrored looking surface. If the quality of the surface is not as described, insert the rod into the sharpener and repeat the cutting and polishing procedure. Remove the rod, inspect the surface quality and if acceptable, place the rod electrode into the original box for storage until ready for use. To prevent contamination of a sharpened rod electrode, do not touch the tip or edge of the tip of the sharpened electrode with the fingers or metallic surfaces or anything but a fresh, clean laboratory grade paper towel. Do not use a rod electrode for analysis if the surface appears to have been damaged. Refer to Section 3.4 of this manual for the procedure to change/rotate the cutter blade.

The electrode sharpener power is on a timing circuit and will turn off after approximately 3 to 5 minutes. The electrode sharpener can also be turned off by placing the MODE switch in the STANDBY position.

1.3 Installing the Disc Electrode

The disc electrode is the most significant contributor to the accuracy and repeatability of the instrument. They are manufactured and then purified to strict specifications to ensure that they do not contain unacceptable levels of trace element contamination for the elements of interest. The

care taken to properly install the disc on the shaft will help to ensure that excitation parameters will be kept as constant as possible, thus resulting in repeatable analytical data.

To install the disc electrode on the shaft, a laboratory grade disposable towel is recommended. See Section 1.7 for a description of the laboratory grade paper towel. Pour out a few disc electrodes onto a clean laboratory grade paper towel. Take a laboratory grade paper towel and double it to be sure that no contamination from the fingers will be absorbed into the disc. If large size paper towels are used, they should be cut with scissors into two inch squares to facilitate easy handling. Place the towel over the disc electrodes, and with the forefinger and thumb, grab one disc from the pile and place the disc on the shaft, Figure 1-4. With firm pressure, push the disc electrode onto the shaft until it comes to rest against the index shoulder of the shaft. If the disc electrode does not offer some resistance to the shaft as it is being inserted, remove and discard this electrode because the inner diameter has not been made to the tolerances specified.

CAUTION: *The disc electrode shaft is designed to be replaced by the operator using a small jeweler's screwdriver. The shaft has right hand threads for tightening it into the commutator. When pushing the disc electrode on the shaft, do not apply counter-clockwise rotation on the disc electrode as this may cause the disc electrode shaft to loosen.*

NOTE: *Loose disc electrodes will produce erroneous results. If the disc electrode is too loose, arcing will occur between the inner diameter of the disc and the outer diameter of the shaft.*

CAUTION: *The disc electrode may be very hot to the touch.*

Use a towel to remove a disc electrode from the shaft after an analysis, and to wipe away any oil which may have spilled over from the burn.

1.4 Installing the Rod Electrode and Setting the Gap

The rod electrode is installed after the disc electrode is already in place. To install the rod electrode, take the rod in the fingers of the right hand and with the left hand apply inward pressure to the round black rod electrode clamp knob, Figure 1-4. This will open the clamp door approximately 3/8 inch (9.5 mm).

Insert the rod electrode into the vertical "v" shaped channel until the sharpened tip can be seen protruding from the bottom of the rod holder and gap setting device. Release the rod clamp knob and the rod electrode will be pinched between the back of the rod clamp knob and the centering "v" channel. Press and then release the rod clamp knob again and the rod electrode will drop by gravity and come to rest on the disc electrode.

Raise the analytical gap setting lever. This action will drive the rod electrode holder and slide mechanism downward along the vertical axis. As the slide mechanism moves downward, the rod electrode remains in the installed position, because there is zero clearance between the disc and rod electrodes. The analytical gap setting lever will reach the end of its travel when it is raised to the full upward position. Return the analytical gap setting lever to the lowered position. As the lever begins to return to the lowered position, the rod electrode holder and slide mechanism begins to raise upward along the vertical axis until it is stopped by the analytical gap adjustment screw. This time the rod electrode, which is clamped in the rod electrode holder and slide mechanism, will travel upward with the slide mechanism. An analytical gap distance of 0.090 inches has now been precisely set.

NOTE: *Care must be taken not to touch the brass block with the tip of the carbon rod in order to avoid false copper readings.*

The analytical gap distance has been set during factory calibration and should not be readjusted




Figure 1-4, Sample Stand and Controls for Routine Operation

during routine operation.

1.5 Installing and Positioning the Sample Holder

The Spectroil M can accommodate several different types of sample holders and a sample holder cover. The following paragraphs describe the procedure to install disposable and reusable sample holders and the sample holder cover.

1.5.1 Disposable Sample Holder

Installing the oil or fuel sample to be analyzed should be the last step in loading the sample stand for analysis. When performing fluid analysis, an important consideration which has an effect on the reproducibility of the analysis is the quantity of the sample introduced into the analytical gap. This parameter is one for which the instrument cannot adjust. Proper level of oil in the sample holder is, therefore, part of any good operator technique. Standards and samples are analyzed either in disposable plastic sample holders or a reusable sample holder. In either case, it is recom-

mended that the sample holder be filled level with the top.

NOTE: An adapter may be required with some of the commercially available disposable sample holders.

With the forefinger and thumb, pick up the sample holder and place it in the slot at the top of the sample table. Push the sample holder towards the back of the sample table until the sample holder comes to a stop. It is now properly positioned in the sample table. Lift the table positioning lever, located on the bottom of the sample table, upward until the table reaches the end of its travel. The bottom of the disc electrode should now be immersed in the sample. The sample is now ready for analysis. Close the door of the sample stand and press the START button located on the operator's control panel or function key 9 (F9) on the keyboard, or the burn icon.

Upon completion of the analysis, open the sample stand door, lower the table and remove the sample holder. Please note that proper oil disposal proce-

dures must be followed as dictated by local regulations and laws.

The sample holder may be HOT to the touch depending on the type of holder and the oil or fuel that was analyzed.

1.5.2 Reusable Sample Holder

The procedure is identical to the above disposable sample holder procedure except that a clean reusable sample holder is used to hold the oil or fuel for the analysis. The sample stand table also has a special cutout and groove to hold the sample holder in place and to align it properly for the analysis, Figure 1-5.

Upon completion of the analysis, open the sample stand door, lower the table and remove the sample holder. Pour the oil into a suitable container for proper disposal and set the sample holder aside for cleaning. Please note that proper oil disposal procedures must be followed as dictated by local regulations and laws. The sample holder should be cleaned with an ultrasonic bath and an environmentally acceptable cleaning solution.

1.5.3 Sample Holder Cover

Some fuel samples and hydraulic oils may catch fire at some point during the analysis. For such samples, a sample holder cover should be used to retard the flame and minimize smoke which will attenuate the signal from the analysis. The cover works only with the reusable sample holder.

The following sample stand preparation sequence should be followed to analyze samples that require the sample holder cover:

1. Install the disc electrode, Section 1.3
2. Install and position the reusable sample holder, Section 1.5.2, Figure 1-5. Raise the sample holder in position and with a disposable pipette, fill the sample holder with the fuel sample. Do not overfill the sample holder.




Figure 1-5, Sample Stand with Reusable Sample Holder in Place

3. Place the cover over the reusable sample holder and disc electrode. Note that the cover only fits in one direction and has a cutout for the disc electrode shaft.
4. Install the rod electrode and set the gap, Section 1.4.

The sample is now ready for analysis, Figure 1-6. Close the door of the sample stand and press the START button located on the operator's control panel or function key 9 (F9) on the keyboard, or the burn icon.

When the analysis is complete, open the sample




Figure 1-6, Sample Stand with Reusable Sample Holder and Cover in Place

stand door, remove the rod electrode, remove the cover, lower the sample table, remove the sample, and remove the disc electrode. The sample stand is now ready for the next analysis.

1.6 Cleaning the Sample Stand

The Spectroil M incorporates the rotating disc arc emission technique for excitation of the fluid sample. This technique produces a fine carbon residue which, when combined with oil droplets, produces an oil coating over the sample stand and door area. If allowed to accumulate, this coating will collect the carbon particles and eventually produce a lower resistance path than the analytical gap. If this occurs, the high voltage will not discharge across the analytical gap, but will discharge along the lower resistance path causing damage to the sample stand components.

To prevent arc-over, it is recommended that the operator perform the simple cleaning procedures outlined below.



WARNING

PROLONGED CONTACT WITH SOME SOLVENTS AND OILS MAY CAUSE CANCER!



WARNING

DO NOT USE ANY CHLORINATED SOLVENTS INTERNALLY OR EXTERNALLY ON THE INSTRUMENT!

CAUTION: All chemicals should be used in accordance with good laboratory practice. Proper ventilation is required when using any solvent. Skin contact and prolonged exposure to fumes produced by any solvent may be hazardous.

1.6.1 Cleaning After Each Burn Cycle

Take the paper towel used to remove the disc electrode from the shaft and clean the shaft, the sam-

ple table, and the sample plate area between the disc electrode shaft and the rod electrode clamp.

1.6.2 Cleaning After Each Operating Shift

After 8 hours of operation, the complete sample stand area must be wiped clean of the oil film buildup created by the burn cycles. If performed routinely, the sample stand can be cleaned simply with paper towels and moderate rubbing. However, if this procedure is performed sporadically or inadequately, an oil dispersant may be required to remove the buildup. A general purpose foam type spray detergent is recommended to dissolve the oil film buildup. A spray detergent is capable of contacting those areas which are hard to reach. Remove all detergent by wiping dry with paper towel.

1.6.3 Cleaning the Quartz Window

The quartz window that protects the lens and fiber optic must be cleaned frequently depending upon the type of fluid being analyzed. In general, this should be done at least every 5 burn cycles. To clean the protective quartz window, take a clean, soft, disposable laboratory tissue and wet one corner of the towel with isopropyl rubbing alcohol or ammonia based window cleaner. With the forefinger, rub the wetted portion of the paper towel along the surface of the window while applying moderate clockwise pressure on the window. This will disperse the oil film. Now take the dry portion of the paper towel and repeat this procedure until no oil can be seen on the tissue paper. A cotton swab can also be used for this purpose. A diluted solution of ammonia and water may be used.

CAUTION: The lens protected by the window does not require cleaning and should only be disassembled by a qualified engineer.

CAUTION: Do not use solvents to clean this window as they may selectively block or attenuate the passage of light necessary to determine the presence and concentration of the elements in the oil samples.

1.6.4 Cleaning Solutions

The Spectroil M is designed to analyze petroleum and synthetic base products. In operation, the handling and actual analysis of these products create spillage and often leave an oily film on the instrument. In general, these spills can be adequately cleaned simply by wiping the surface with paper towels. There are occasions, however, where the petroleum/synthetic product may require a detergent to dissolve the petroleum base. For these occasions and for routine cleaning, a general purpose spray and wipe detergent is recommended for internal and external instrument components.

1.7 Paper Tissue for Operating and Paper Towels for Cleaning

Disposable paper tissues and towels are recommended for use in the daily operation of the Spectroil M. The type of paper tissue used to handle the disc electrodes is very important. Most household tissue paper is treated with certain elements to make it softer or more absorbent. If used to handle the disc electrodes, these elements will contaminate the electrodes and produce erratic results, especially for silicon. Therefore, it is recommended that a laboratory grade paper tissue be used for this operation.

Paper towels are useful to clean the sample stand components and wipe spills which occur during routine operation. The type of paper towel used for this function is not critical. Typical household towels or C-fold janitorial towels work best for this function because of their absorbent characteristics.

1.8 Waste Oil Disposal Container

A waste oil container for oil analysis applications is required to properly dispose of the remaining oil sample after the analysis cycle. It is recommended that this waste oil container be in the form of a rectangular pan approximately 6 inches long, 4 inches wide, and 1 inch deep, with a

screened cover to permit the remaining oil to drain through the screen. If a drain tube is installed on the bottom of the waste oil container, the waste oil container can continuously empty into a large capacity reservoir for proper disposal. Good laboratory procedures should be exercised in the disposition of all waste oils.

2.0 DAILY OPERATION

This section details those procedures that will be routinely used in the day-to-day operation of the Spectroil M. The operator must be familiar with the general operating requirements described in Section 1.0. A flow chart of the normal daily routines is shown on the next page in Figure 2-1. For convenience, the parentheses after each step in the chart refer to the corresponding sections in this chapter. The various procedures are explained in brief, easy to follow step-by-step instructions.

2.1 Daily Routine Prior to Use

1. Place the MODE switch on the operator's control panel to the OPERATE position, Figure 2-2. Power will be applied to the electrode sharpener and exhaust fans when MODE switch is in the OPERATE position.
2. Verify positive action from the sample stand exhaust system. With the sample stand door open, hold a piece of tissue paper up to the exhaust filter. It should be sucked up and held in place against the filter. Remove the



Figure 2-2, Control Panel




Figure 2-1, Daily Operating Procedure Flow Chart

tissue.

3. Turn printer ON and check to see that sufficient paper is available. If the printer has an ON LINE light, it should be illuminated.
4. Ensure that an ample supply of sample holders, sharpened electrodes and discs are on hand.
5. Select standards for daily use and shake vigorously for at least 30 seconds.
6. Have an oil waste container on hand (Section 1.8).
7. Have cotton swabs, contaminant free tissue paper and paper towels on hand (Section 1.7).

2.2 Warm-Up Procedure

If the Spectroil M has been idle for several hours, it may be necessary to conduct a series of burns to introduce light into the optics and to allow the electronics to become warm. This warm-up exercise can be conducted with any oil sample or standard and can use electrodes which have been burned before. It is recommended that at least three warm-up burns be conducted.

1. Analyze or “burn” three or four samples (do not burn the same sample more than twice to prevent sample ignition) in accordance with the instructions given in Section 2.3 Routine Sample Analysis.

NOTE: For the warm-up cycle only, the same disc and rod electrodes can be utilized for up to four consecutive burns but the electrodes have to be re-gapped after each one.

2. The results produced by the warm-up burns are of no use. Press function key 6 (F6) AVERAGE . This prepares the screen for the next analysis at which time the three warm-up burns will be cleared.

2.3 Routine Sample Analysis

This paragraph gives the steps to follow to analyze or “burn” any type of sample, whether it is a used oil sample, an oil standard, coolant, or a fuel sample. Refer to the referenced sections for details. The various parts referred to are shown with labels in Figure 1-4.

NOTE: When a new lot of disc electrodes is started, either from a new manufacturer or a different lot from the same manufacturer, the disc electrode offset procedure of Section 3.3 must be performed.

1. The video monitor should display the Analysis Program screen, Figure 2-3. If a screen saver is in use, the Analysis Program screen will not be displayed. Press any key on the keyboard to terminate the screen saver and re-display the Analysis Program screen.
2. Install a carbon disc on the disc shaft using a clean laboratory grade tissue to avoid contact with fingers (Section 1.3).
3. Press inward on the black plastic knob of the spring loaded rod electrode clamp to open the jaws of the clamp. Insert a graphite rod electrode until the tip of the carbon rod is in contact with the disc electrode, then release the knob to secure the electrode in the clamp.




Figure 2-3, Sample Analysis Screen

4. Set the analytical gap mechanism by raising and then lowering the analytical gap setting lever. This will set a gap distance of 0.090 inches between the disc and rod electrodes (Section 1.4).
5. Fill a sample holder with sample to be analyzed. Be sure to always fill sample holders to the rim (Section 1.5).

NOTE: When light fuel samples are analyzed, reusable sample holders and sample holder covers must be used to prevent sample ignition.

6. Place the filled sample holder on the table and slide it back to the end of the groove on the table (Section 1.5).
7. Raise the sample table using the sample table positioning lever (sample fluid will contact bottom of disc) (Section 1.5). Press function key 3 (F3) if sample ID's are to be entered, or see Section 2.8 to set-up sample ID's.
8. Close the sample stand door and press the START button or function key 9 (F9) START.
9. After the burn is complete, results will appear on the video screen.
10. After the burn is complete, open the sample stand door and remove the rod electrode. Set it aside for subsequent re-sharpening before it is used again.
11. Lower the sample table, remove and discard the sample holder.

CAUTION: *The disc electrode will be hot to the touch.*

12. Using a paper towel to protect fingers from the hot disc, remove and discard disc electrode.

13. Using a tissue or paper towel, wipe excess, spilled or splattered sample fluid from sample table and disc electrode shaft.

NOTE: The quartz protective window should be cleaned at least every 5 burns.

The Spectroil M is now ready to analyze another sample by again following the above 13 steps.

2.4 Daily Standardization Check

The standardization check is performed to verify that the instrument has remained in calibration. It is a quick method of verifying that the instrument can give accurate results without conducting a full standardization.

This procedure requires that the operator analyze at least three different levels of calibration standard. One of the standards should be a base oil or 0 ppm standard, the next standard should be at the high end of the concentration range expected in the unknown samples, and the third standard should be some concentration between the 0 ppm and the high standard. For example, if the samples to be analyzed are used oils which normally have iron concentrations as high as 100 ppm and silver concentrations as low as 1.0 ppm or less, the recommended standards for the daily standardization check should be the base oil (0 ppm) and the 100 ppm standard. These two standards will cover the complete calibrated range from 0 ppm to 100 ppm for all elements. The third standard may be 10 or 30 ppm.

For military aircraft applications, the 0, 5, 10 and 30 ppm standards generally cover the complete range of expected concentrations. If it is a fuel sample with very low contamination, the recommended standards are the base oil and the 10 ppm standard oil. The following steps should be conducted only after the warm-up procedure has been completed and the window has been cleaned in accordance with Section 4.1.6.3.

1. Prepare sample stand in accordance with Daily Routine Prior to Use, Section 2.1.
2. Make 3-4 warm-up burns of D12-100 PPM standard.
3. Choose three standards in the expected range of the samples.
4. Clean the quartz window.
5. Make one burn of a standard selected in step 3 above.
6. Compare the results of this burn with Table 2-1. If all the elements are within the acceptable range, proceed to Step 8, 9, or 10. If the results are not within the range, proceed to Step 7.
7. Make a two more (total of three) burns of the standard selected in step 5 above and press function key (F6). If all elements are within the range, proceed to step 8, 9, or 10. If not, perform complete Standardization Procedure in accordance with the Daily Operating Procedure Flow Chart.
8. Choose the second standard from step 3 and repeat steps 4- 7.
9. Choose the third standard from step 3 and repeat steps 4-7.
10. Daily Standardization Check is now complete.

NOTE: Table 2-1 provides a recommended range that all elements should fall between during a daily standardization check. The ranges are narrow because they are based on one or three analyses and should not be confused with the actual accuracy and repeatability specifications for the spectrometer given in Tables 4-1 and 4-3 which are based on ten analyses. If Table 2-1 ranges can be met, then it is assumed that by default, Tables 4-1 and 4-2 will be satisfied.

Table 2-1, Acceptable Range Indices for Daily Standardization Check

Concentration	Min.	Max.
0	0.0	1.0*
5	3.8	6.2
10	8.5	11.5
30	27.0	33.0
50	45.0	55.0
100	90.0	110.0
300	255	345

* This range applies to all elements except Ag, Al, Mg and Cu. The range for these elements is 0 to 0.5.

2.5 Complete Standardization

Complete standardization is a procedure performed to place the calibration of the instrument as close to the standard values as the instrument originally produced during factory calibration. This procedure involves burning oil standards at predetermined points along the calibration curve. After these standards are analyzed, the computer software will determine mathematical factors to correct for any change in the calibration. Complete standardization is performed under the following conditions:

- When the instrument has been relocated to another site for operation. This is generally performed after the optical profile procedure has been completed.
- When results from the daily standardization check fall outside of acceptable limits for operation.
- Prior to the analysis of JOAP monthly correlation samples.
- After optical profiling procedure has been performed.

A complete standardization is performed by burning two or more calibration standards that have been pre-selected during factory calibration of the instrument. The concentration levels for complete standardization have been selected based on the application and typical operating range for the elements of interest. In general, all

elements are standardized at 0 ppm to determine the background level, all wear metal and contaminant elements are standardized at 100 ppm and additive elements on commercial instruments are standardized at 900 ppm. These concentration levels are programmed into the computer and are displayed at the appropriate time in the following procedure.

From the Analysis Program screen, press function key 7 (F7), choose the Standardization icon, or select the Operations/Standardize pull down menu. The software will automatically clear all previous measurements from the video display. A dialog with the name of the first calibration standard the instrument will expect to measure will appear. Refer to Figure, 2-4.

Three options exist when this dialog appears. The first is to select the OK button which confirms that the operator will begin to measure the 0 PPM standard. The second option is to select the SKIP button indicating the operator does not intend or need to measure the 0 PPM standard and wants to increment to the next standard in the standardization process. The last option is to press CANCEL and this action will terminate the standardization routine completely and return the software to the Analysis Program screen.

In most cases, the OK button will be selected and the dialog will disappear. Centered just below the tool bar will be the name of the standard, 0 PPM in a red banner. This banner will remain there




Figure 2-4, First Standardization Sample Dialog

until an average is made and the next standard will appear. Those elements to be standardized at this concentration level will have their values appear and those elements not standardized will have no values appear.

NOTE: Reference channels are not standardized and therefore will not appear highlighted.

1. Select the 0 ppm standard and fill five sample holders. Take care to always fill the sample holder to the rim.
2. With an optical lens cleaning solution (not containing silicon), isopropyl alcohol, or a window cleaning solution with ammonia, clean the quartz window attached to the Fiber Optic Lens Holder in accordance with Section 1.6.3.
3. Following procedures set forth in Section 2.3, burn all five samples of the 0 ppm standard.
4. On completion of Step 4, look at the readings on the video screen. If one of the five burns does not appear to represent the other four, it may be rejected.

NOTE: The decision to accept or reject burns during this procedure is at the discretion of the operator. Quite often an operator will know the cause of a rejectable burn and therefore reject it almost automatically. Rejectable burns can be caused by inconsistencies in consumables such as excess variation in the specific density (hardness) of the disc electrode, a loose fitting disc electrode, a poorly sharpened rod electrode, an under- or over-filled sample holder or an analytical gapping error.

To help as a guideline in making a determination whether to reject a burn or not, it is recommended to follow the 80% rule. This rule states that a burn qualifies to be rejected if 80% or more of the elements exhibit the same symptom. For example, a rejectable burn is one which is obviously too high or too low when compared to the other four

burns. Take for example five burns of 100 ppm for the element Fe (100, 102, 115, 107 and 99), Ag (107, 110, 124, 108 and 106), Al (98, 95, 96, 99, and 93), etc. In this example, the third burn is suspect of being a high burn. To determine if the third burn qualifies to be rejected, count the elements in which the third burn was the highest of the five measurements. In the abbreviated example above, Fe and Ag meet this criteria but Al does not. Total the elements that exhibit this condition and if 80% (8 out of 10) meet this criteria, this burn qualifies to be rejected.

To reject one of the measurements from the video display, move the selection pointer over any portion of the measurement and left click the mouse one time. This will highlight the burn in a black background. Pressing the DELETE key on the keyboard one time will remove this measurement off of the screen. In the event that the wrong measurement has been highlighted, position the pointer over the measurement again and left click the mouse a second time. This will remove the highlight, then select the proper measurement to reject and press the DELETE key.

If a group of sequential measurements are to be rejected, for example measurements 3 through 7, place the pointer over the first measurement (#3) and left click the mouse one time to highlight the measurement. Then move the pointer over the last measurement (#7) and hold the SHIFT key down and left click the mouse one time. All measurements from number 3 through number 7 will be highlighted. Press the DELETE key to remove all five measurements.

5. After making as many measurements as necessary to obtain a good average, press function key 6 (F6), click the average icon, or select Operations/Average from the pull down menu. The average of each element will be calculated and displayed below each element's column of measurements. To make a printout of the measurements and their average, press the print icon. Automatically, the next dialog will appear provid-

ing instructions to make measurement of the next standard. The name of the standard will be different from instrument to instrument depending on application (M, M/N, M/F, M/C) and customer specifications. In general, the next standard will be either a 10/30 PPM or a 100 PPM concentration as shown in Figure 2-5.

6. Select the next standard (i.e., 100 PPM) and fill five sample holders. Take care to always fill the sample holder to the rim.
7. With an optical lens cleaning solution (not containing silicon), isopropyl alcohol, or a window cleaning solution with ammonia, clean the quartz window attached to the Fiber Optic Lens Holder in accordance with Section 1.6.3.
8. Following procedures set forth in Section 2.3, burn all five samples of the next standard (i.e., 100 PPM).
9. On completion of Step 8, look at the readings on the video screen. If one of the five burns does not appear to represent the other four, it may be rejected. Refer to the note regarding the 80% rule to determine if a burn qualifies to be rejected.
10. After making as many measurements as necessary to obtain a good average, press function key 6 (F6), click the average icon,




Figure 2-5, Second Standardization Sample (Standard 2)

or select Operations/Average from the pull down menu. The average of each element will be calculated and displayed below each element's column of measurements. To make a printout of the measurements and their average, press the print icon. Automatically, the next dialog will appear providing instructions to make measurement of the next standard. The name of the standard will be different from instrument to instrument depending on application (M, M/N, M/F, M/C) and customer specifications. In general, the next standard will be either a MA 900 PPM or an AM SPECIAL1000 PPM concentration as shown in Figure 2-6.

11. Select the next standard (i.e., MA 900 PPM OR AM SPECIAL 1000 PPM) and fill five sample holders. Take care to always fill the sample holder to the rim.
12. With an optical lens cleaning solution (not containing silicon), isopropyl alcohol, or a window cleaning solution with ammonia, clean the quartz window attached to the Fiber Optic Lens Holder in accordance with Section 1.6.3.
13. Following procedures set forth in Section 2.3, burn all five samples of the next standard (i.e., MA 900 PPM).
14. On completion of Step 13, look at the readings on the video screen. If one of the




Figure 2-6, Third Standardization Sample (Standard 3)

five burns does not appear to represent the other four, it may be rejected. Refer to the note regarding the 80% rule to determine if a burn qualifies to be rejected.

15. After making as many measurements as necessary to obtain a good average, press function key 6 (F6), click the average icon, or select Operations/Average from the pull down menu. The average of each element will be calculated and displayed below each element's column of measurements. If there are no additional standards to be measured as part of the standardization routine, next dialog to appear will indicate standardization is complete and inquire if the average and burns for the last standard measured should be printed. Refer to Figure 2-7 below.
16. Click on the Yes button if a printed copy of the analyses is desired. After clicking on the Yes or the No button, the Standardization Values Screen, Figure 2-8 is displayed. This table can also be displayed at any time by selecting Program/Standardization Samples/Values.

NOTE: It is strongly recommended that accurate records of the complete standardization data be kept for future reference. For this reason we highly recommend that printouts of all daily standardizations be printed and kept on file. This data reflects the current Spectroil M standardization.




Figure 2-7, Standardization Complete, Dialog

The table in Figure 2-8 compares the expected intensities (the intensities generated during factory calibration) to the obtained intensities (the intensities obtained from the most recent complete standardization). From the relationship of the expected intensities compared to the obtained intensities of the low standard and the high standard, the standardization factor is derived.

Another table is displayed after you click OK or Cancel on the Standardization Values Screen, Figure 2-9. This second table is called the Standardization Factors Table. This table can also be displayed on the screen by selecting Program/Standardization Samples/Standardization Factors.

This table displays a factor for each element. At

Element	Low Sample	Expected	Obtained	High Sample	Expected	Obtained
1. Ag	0 PPM	0	43	100 PPM	664000	629422
2. Ag2	0 PPM	9952	9988	100 PPM	18453	18196
3. Al	0 PPM	1	8	100 PPM	167000	152637
4. B	0 PPM	1	51	100 PPM	143000	142390
5. Ba2	100 PPM	44835	40321	MA 900 PPM	234000	214952
6. Ba	0 PPM	1	590	100 PPM	8833000	7738775
7. Ca2	0 PPM	230000	233072	100 PPM	414090	400810
8. Ca	0 PPM	3988	2457	100 PPM	11516000	11510745
9. Cd	0 PPM	31	68	100 PPM	60231	58220
10. Cr	0 PPM	280	69	100 PPM	1088000	913189
11. Cu2	0 PPM	1	111	100 PPM	6064	7448
12. Fe	0 PPM	51	56	100 PPM	57325	52135
13. H						
14. Mg	0 PPM	126	191	100 PPM	306000	283868
15. Mn	0 PPM	12	46	100 PPM	170000	143990
16. Mo	0 PPM	6	22	100 PPM	32757	31231

Figure 2-8, Standardization Values Table (Example)

Element	Low Sample	High Sample	Factor	Offset
1. Ag	0 PPM	100 PPM	1.055	53
2. Ag2	0 PPM	100 PPM	1.084	-1279
3. Al	0 PPM	100 PPM	1.034	-8
4. B	0 PPM	100 PPM	1.004	52
5. Ba	0 PPM	100 PPM	1.115	659
6. Ba2	100 PPM	MA 900 PPM	1.083	1158
7. Ca	0 PPM	100 PPM	1.000	1540
8. Ca2	0 PPM	100 PPM	1.090	-23062
9. Ca3	100 PPM	MA 900 PPM	0.914	24299
10. Cd	0 PPM	100 PPM	1.015	98
11. Cr	0 PPM	100 PPM	1.191	198
12. Cr2	0 PPM	100 PPM	0.929	38
13. Cu	0 PPM	100 PPM	1.074	151
14. Cu2	0 PPM	100 PPM	1.067	119
15. Fe	0 PPM	100 PPM	1.100	-11
16. Mg	0 PPM	100 PPM	1.078	-80

Figure 2-9, Standardization Factors Table (Example)

best, the factor would be exactly 1.000. In practice, differences in electrode grades, standards, and instrument variables will cause the intensities achieved to result in factors which are either slightly above or below 1.000. These factors should remain somewhere between 0.5 and 5.0. If a factor exceeds these tolerances, it is not always an indication of an error or pending problem. If such a case should occur, consult Spectro Incorporated Field Service for analysis and explanation.

18. A daily standardization check in accordance with Section 2.4 should be carried out to verify calibration.
19. Standardization is now complete and it is possible to burn routine used oil samples.

2.6 Daily Routine Prior to Securing

1. Turn the MODE switch on the control panel to STANDBY.
2. Turn the printer power switch OFF.
3. Remove disc and rod electrodes.
4. Clean disc electrode shaft with a paper towel.
5. Clean and wipe the entire sample stand area.
6. Clean the quartz protective lens using a clean soft disposable laboratory tissue.
7. Wipe all oil standard bottles clean.
8. Check supply of standards (don't run out).
9. Clean and wipe used oil container.
10. Clean the working area.
11. Sharpen all rod electrodes and store them so they are protected from inadvertent contamination.

2.7 Optical Profiling

The Spectroil M optical system is shock mounted in a light-sealed and environmentally protected temperature stabilized enclosure. Consequently, the optics do not need to be profiled frequently. However, detection limit and repeatability suffer when the optics are off profile. Unfortunately, there is no one rule which ensures that the optics are on profile. The following guidelines are presented to indicate when profiling should be done:

- At least once every month prior to analysis of JOAP monthly correlation samples.
- After the instrument has been transported to a new location.
- Whenever the instrument has been subjected to temperature variations greater than 15° F (10° C).

If one of these apply, it is also reasonable to perform a standardization as detailed in Section 2.5.

Follow the next steps as detailed in this procedure to determine the optimum optical peak profile position for operation.

1. Prepare the sample stand with new electrodes in accordance with the appropriate paragraphs of the General Operating Requirements section of this manual.

NOTE: it is important to only a 100 PPM standard for the profiling procedure (D12-100 for military instruments and 100 PPM standard for commercial instruments)..

1. To begin the optical profile procedure press function key 4 (F4), left click the profile icon, or select Operations/Profile from the pull down menu options.
2. The screen shown in Figure 2-10 appears and calls for the first analysis of the profiling standard..
3. Fill a sample holder with the profiling stan-




Figure 2-10, Optical Profile Screen

dard and following procedures set forth in Section 2.3, to burn an oil sample.

4. At the completion of the burn, a screen similar to Figure 2-11 will appear with the profile log. The profile log shows the current profile for each chip, the previous profile and the difference between the two. The status for each chip should be OK, if not repeat the profile procedure one or two more times until the status is OK for all chips. If this condition cannot be achieved after three attempts, contact Spectro for assistance.
5. Click on OK to return to the main analysis screen. The profiling procedure is now complete.

2.8 Sample Identification (I.D.)

A sample identification (I.D.) can be added to each analysis for identification purposes. The

Chip ID	Offset	Status	Last	Status	Diff	Status
1	0.344	IDE HWB TOO LARGE	0.344	IDE HWB TOO LARGE	0.000	OK
2	-5.564	IDE HWB TOO SMALL	-5.564	OK	0.000	OK
3	-0.627	IDE NO PEAK	-0.627	IDE NO PEAK	0.000	OK
4	-2.355	OK	-2.424	OK	0.039	OK
5	-0.141	IDE NO PEAK	-0.141	IDE NO PEAK	0.000	OK
6	-9.875	OK	-9.894	OK	0.019	OK
7	-2.551	OK	-2.530	OK	0.079	OK
8	-0.077	IDE HWB TOO SMALL	-0.077	IDE NO PEAK	0.000	OK
9	-3.300	OK	-3.309	OK	0.009	OK
10	0.442	IDE NO PEAK	0.442	IDE NO PEAK	0.000	OK
11	0.544	IDE NO PEAK	0.544	IDE HWB TOO SMALL	0.000	OK
12	0.100	IDE UNDEF CHIP	0.100	IDE UNDEF CHIP	0.000	OK
13	-0.591	IDE HWB TOO SMALL	-0.591	OK	0.000	OK
14	1.375	OK	1.382	OK	-0.007	OK
15	-6.412	OK	-6.400	OK	-0.012	OK
16	5.191	OK	5.233	OK	-0.102	OK

Figure 2-11, Optical Profile Log

ample I.D. can be configured by the user to include a variety of information about the analyzed sample. This section describes the setup process to format the sample I.D. and how to enter it in routine operation.

2.8.1 Setup of the Sample I.D.

To configure the format for the sample I.D., select "Sample I.D." from the System pull-down menu. This menu option produces a dialog, Figure 2-12, that can be configured by the user to meet any combination of alpha-numeric characters for global sample identification.

A sample can be identified by up to six field segments. This dialog permits the operator to choose how many field segments will be used for the sample identification and name each of these field segments.

Each field segment can also be auto-incremented, which means that after the first sample number is entered, and if all numbers that follow are in numeric order, they can be automatically filled in ascending order thus saving time.

The sequence of how the fields will appear can be determined or altered and the size of each field can be customized up to a maximum of 40 characters total. The last column of the sample I.D. dialog is for V4.2 protocol. This protocol is capable of storing up to two segments of the sample ID in the file. More than three segments are not permitted. That portion of the sample identifi-

cation that has V4.2 assigned as segment 1 or segment 2 will be transmitted or stored under the V4.2 protocol.

2.8.2 Using Sample ID's

When the Spectroil M is ready to analyze a sample, the Sample I.D. can be entered by pressing function key 3 (F3), or by clicking on the Sample I.D. icon. This menu option will produce a dialog to permit the entry of one single ID, Figure 2-13, or provide the capability to pre enter multiple sample ID's, Figure 2-14, using the MULTIPLE button. Both dialogs are configured at the system level by the System/Sample ID option.

The Multiple Sample ID Entry dialog allows up to 50 sample identifications to be pre loaded to facilitate rapid sample throughput. All sample ID fields are configured at the system level through the System/Sample ID menu option. Along the bottom of the Multiple ID dialog are buttons to expedite the entry of sample numbers. The Copy button will copy the contents of one field and permit it to be copied into another field of equal or greater field size using the Paste button. Copy All will copy one sample and insert it into all remain-




Figure 2-13, Single Sample ID Entry Menu




Figure 2-12, System/Sample ID Setup Menu




Figure 2-14, Multiple Sample ID Entry Menu

ing empty fields in that column. If that column is set for auto increment, the sample number will increase one value per row. Insert, Delete, and Clear are self-explanatory.

Click load to proceed. The first time that this option is enabled, an input file layout screen, Figure 2-15 appears and must be filled in with the user layout preferences. Click OK when complete and the a screen enabling the user to select the sample ID files appears, Figure 2-16.

2.9 Sulfur Analysis

An optional sulfur analysis capability is available for the spectroil M/C and M/N. The optional capability consists of a separate sulfur optic attached to the right hand side of the sample stand and some minor modifications to the sample stand. The sulfur spectral line is in a region of the spec-

trum where any signal resulting from the excitation process is absorbed by air. For this reason, when sulfur analysis is required, the sample stand area and the optic must be purged with an inert gas so that the CCD detector is able to see the sulfur signal generated at the analytical gap. To solve this problem, the Spectroil M uses nitrogen as a gas and a separate optic just for the sulfur capability.

2.9.1 Set-up for Sulfur Analysis

Follow this procedure to prepare the sulfur optic for routine analysis.

1. Attach a low pressure regulator onto a Nitrogen bottle or other nitrogen source. The output of the regulator should have a 1/4 inch barbed quick disconnect, or ferrule hose fitting
2. Attach the 1/4 inch ID (inner diameter) hose that comes out of the back of the sulfur optic to the 1/4 inch barbed quick disconnect, or ferrule hose fitting on the regulator. Verify that all hose fittings are secure and no pressure leaks are present. Teflon tape should be used on all threaded hose fittings.

CAUTION: Pressure leaks will cause the introduction of fresh air to the purge path and will result in low readings for sulfur.

Sulfur is analyzed separately from the other elements, therefore the sample stand must be set-up so that the excitation process is viewed by the sulfur optic and not the routine optic. The excitation process can be selectively viewed to the left by way of the sulfur light pipe, or to the right for the routine analysis by way of the lens assembly.

3. To activate the view to the sulfur optic, the "L" bracket must be removed from the sulfur optic light pipe. Unscrew the round knurled nut that holds the "L" bracket in place and remove it from in front of the sulfur light pipe, Figure 2-17




Figure 2-15, Sample ID File Layout Preference Screen




Figure 2-16, Sample ID File Selection Screen




Figure 2-17, "L" Bracket Installed Covering Sulfur Optic Light Pipe

4. Relocate the "L" bracket over the lens assembly of the routine optic located to the left of the analytical gap, Figure 2-18. Use the round knurled nut to fasten it in the new location. This will help to keep the lens assembly clean during the sulfur analysis.

NOTE: Some instruments may be delivered with two "L" brackets installed. When that is the case, loosen the round knurled nut, flip up the "L" bracket for the sulfur light pipe and tighten the round knurled nut. Verify that the other "L" bracket over the lens assembly is in the "down" position.

5. Open the nitrogen bottle to allow the nitrogen to pass through the regulator, adjust the




Figure 2-18, "L" Bracket Installed Covering Lens Assembly

regulator to approximately 25 PSI.

6. Turn the main Nitrogen valve on the sulfur optic control panel to "ON", Figure 2-19.
7. If the Spectroil M is OFF, turn it ON by raising the circuit breaker (CB1) on the right side of the instrument. The spectrometer will power up and the OilM Windows software program will load.
8. Depress the "OPTIC PURGE" button, Figure 2-19, on the sulfur optic control panel for 5 seconds, and release for 5 seconds. Repeat this process 5 times to flush the optical system.
9. To load the Sulfur program select **File** from the pull down menu options, then **Open**. A list of available programs will be displayed. Select the radio button for the Main Sulfur program and press **OK**. The sulfur optic is now ready for operation.

WARNING: Irritating and toxic hydrogen sulfide gas may be found in confined vapor spaces when samples with large percentages of sulfur are analyzed. Greater than 15 - 20 ppm continuous exposure (in air, not to be confused with the liquid calibration standard) can cause mucous



Figure 2-19, Sulfur Optic Control Panel

membrane and respiratory tract irritation. 50 - 500 ppm can cause headache, nausea, and dizziness, loss of reasoning and balance, difficulty in breathing, fluid in the lungs, and possible loss of consciousness. Greater than 500 ppm can cause rapid or immediate unconsciousness due to respiratory paralysis and death by suffocation unless the victim is removed from exposure and successfully resuscitated.

The "rotten egg" odor of hydrogen sulfide is not a reliable indicator for warning of exposure, since olfactory fatigue (loss of smell) readily occurs, especially at concentrations above 50 ppm. At high concentrations, the victim may not even recognize the odor before becoming unconscious.

2.9.2 Profiling the sulfur Optic

The sulfur optic is profiled at the same time as the standard optic. Profiling should be performed if the instrument has been relocated and or shut off for an extended period of time.

The Nitrogen purge should be turned OFF at the sulfur optic control panel during profiling. To profile the optic, follow the instructions in Section 2.7 in this manual.

Note: During the profiling procedure, both optics must see the analytical gap. Make sure the "L" brackets are removed from the sulfur light pipe and the lens assembly. If two "L" brackets are installed, they should be in the up position.

2.9.3 Standardization of the Sulfur Optic

The software must be in the sulfur program to perform this function. If not, see step 9 in Section 2.9.1.

1. From the main analysis screen, press function key 7 (F7), choose the Standardization icon or select the Operations/Standardize pull down menu. The software will automatically clear all previous measurements from the video display. A dialog box with the name of the first calibration standard the instrument will expect to measure will ap-

pear. A blank oil is used on this step and the dialog box may call for "Base Oil", "Blank Oil", or 0 PPM Oil.

2. Following the procedure set forth in Section 2.3, burn three samples of Base/Blank oil and press function key 6 (F6), click the Average icon or select Operations/Average from the pull down menu to average them. To make a printout of the measurements and their average, press the print icon. Automatically, the next dialog box will appear providing instructions to make measurements of the next standard. The concentration of this standard is 1,000 PPM. The name in the dialog box will vary depending on the calibration of the Spectroil M.
3. Following the procedure set forth in Section 2.3, burn the three samples of the 1,000 PPM sulfur standard and press function key 6 (F6), click the Average icon or select Operations/Average from the pull down menu to average them.
4. The next dialog box that appears will indicate standardization is complete and inquire if the average and burns for the last standard measured should be printed, select YES.
5. After confirmation that the burn data and average are to be printed, the last dialog box to appear will inquire if the Factors and Offsets should be printed, again select YES. The Standardization procedure is now complete. Review the factor for sulfur; it should be close to one, +/- 20 %.

2.9.3 Verification

It is good practice to verify calibration after the Spectroil M has been standardized.

1. In PPM mode, run five samples of each calibration standard in the set. Print all burns, averages and statistics. Remember the readings are displayed in percentages.

2. Depending on the type of standard analyzed, the averages should be better than +/- 15%. If not, repeat the standardization procedure, Section 2.9.3 and if it does, proceed with step 3.
3. Routine samples are analyzed for sulfur content per the procedure in Section 2.3, Routine Sample Analysis.

2.9.4 Clean Up and Final Checks

This section provides the steps that must be performed when the sulfur analysis procedure is complete and the Spectroil M/F-LD will be used for the analysis of routine furl samples.

1. Turn the Main Nitrogen valve "OFF".
2. Remove the interchangeable "L" bracket from the lens assembly and place over the end of the sulfur light pipe in the sample stand. This will prevent any oil from traveling down the pipe and contaminating the sulfur optic lens assembly. If two "L" brackets are installed, flip down the sulfur light pipe bracket and flip up the lens assembly bracket.
3. Place the Mode/Operate switch on the Spectroil M/F-LD control panel to **STANDBY**.
4. Thoroughly clean the sample stand area, especially the quartz window and gap sensors. They will be extremely dirty since the nitrogen was blowing the sample flame towards the left side of the sample stand during the sulfur analysis.

2.10 Coolant Analysis

Spectro Incorporated has developed a technique whereby a Spectroil M/C-W or Spectroil M/N-W as applied to used oil analysis can be modified and calibrated to also effectively analyze engine coolant samples; a mixture of approximately

50% water and 50% glycol. The added capability provides the laboratory with a supplementary tool to increase its capabilities and effectiveness. The technique has been shown to be effective for the analysis of wear metals, contamination and supplemental coolant additives in ethylene and propylene glycol.

To be effective, a used coolant analysis program should determine both the coolant condition and the presence of any contaminants or debris. The coolant fluid can be used as a diagnostic medium as the coolant carries not only heat away from the engine parts but also carries fine debris from the interior surfaces of the cooling system. Analysis of the wear debris can provide important information about the condition of the internal parts of the cooling system.

The elements routinely detected and quantified by the Spectroil for coolant analysis are shown in Table 2-2.

2.10.1 Coolant Program Standardization

In order to prepare the Spectroil M/C-W or Spectroil M/N-W for the analysis of coolants, the software must be in the "COOLANT" program.

1. To load the COOLANT program select **File** from the pull down menu options, then **Open**. A list of available programs will be displayed, Figure 2-20. Select the radio button for the GLYCOL (MAIN COOLANT PROGRAM) and press **OK**. The program is now loaded and ready for the standardization procedure.

Table 2-2, Typical Elements Detected in Coolant Samples

<i>Wear Metals</i>	<i>Contaminants</i>	<i>Supplemental Additives</i>
Iron	Silicon	Potassium
Zinc	Magnesium	Silicon
Lead	Calcium	Boron
Copper		Sodium
Aluminum		Molybdenum
Magnesium		Phosphorus




Figure 2-20, Select Program Screen

2. From the main analysis screen, press function key 7 (F7), choose the Standardization icon or select the Operations/Standardize pull down menu. The software will automatically clear all previous measurements from the video display. A dialog box with the name of the first calibration standard the instrument will expect to measure will appear. A blank water sample is used or this step and the dialog box will call for "D.I. Water, Figure 2-21.
3. Following the procedure set forth in Section 2.3, burn three samples of D.I. Water and press function key 6 (F6), click the Average icon or select Operations/Average from the pull down menu to average them. To make a printout of the measurements and their average, press the print icon. Automatically, the next dialog box will appear providing instructions to make measurements of the Glycol standard, Figure 2-22.




Figure 2-21, First Coolant Program Standardization Point




Figure 2-22, Second Coolant Program Standardization Point

The Glycol standard is a multi level standard and the concentrations of the elements vary as stated on its label. A typical makeup of the Glycol standard is shown in Table 2-3.

4. Following the procedure set forth in Section 2.3, burn the three samples of the Glycol standard and press function key 6 (F6), click the Average icon or select Operations/Average from the pull down menu to average them.
4. The next dialog box that appears will indicate standardization is complete and inquire if the average and burns for the last standard measured should be printed, select YES, Figure 2-23.
5. After confirmation that the burn data and average are to be printed, the last dialog box to appear will inquire if the Factors and Offsets should be printed, again select YES. The Standardization procedure is now complete. Review the factors, they should be close to 1, +/- 20 %.

Table 2-3, Typical Elements in Glycol Standard

<i>Concentration</i>	<i>Elements</i>
50 PPM	Al, Ca, Cu, Fe, Mg, Pb, Zn
500 PPM	Mo, Si
1,000 PPM	B, Na(High), K (if installed)
2,500 PPM	P




Figure 2-23, Standardization Complete, Dialog

2.10.2 Verification

It is good practice to verify calibration after the Spectroil M has been standardized.

1. In PPM mode, run five samples of Glycol standard. Print all burns, averages and statistics.
2. Depending on the type of standard analyzed, the averages should be better than +/- 15% for the concentrations listed on the Glycol standard bottle. If not, repeat the standardization procedure, Section 2.10.1.

2.10.3 Routine Coolant Sample Analysis

Coolant samples are analyzed on the Spectroil M/C-W or Spectroil M/N-W as described in Section 2.3, Routine Sample Analysis.

2.10.4 Grease Analysis

Grease analysis with the Spectroil M is similar to analyzing used oil. However, it does present two hurdles that must be considered and overcome. The first concern is obtaining a representative sample of wear metals and/or contamination in the grease, and the second is to physically get the sample into the analytical gap of the spectrometer for analysis.

Obtaining a representative sample is important and tricky for greases. Depending upon how the sample is taken, it may be from the wear track and it will have a high concentration of wear metals, or it may come from a dead space in the mechanical system and the readings will not be much different from a new, unused grease sample. Bear in

mind that grease does not mix in a mechanical system the same way oil does.

One opportunity for sampling may be when a bearing or other grease lubricated component is being charged with fresh grease. It is recommended to collect all of the used grease being forced out of the component by new grease. Then stir the used grease to homogenize it. Now take a sample of the homogenized used grease for further analysis. This procedure is intended to distribute whatever wear and contaminant particles are present as evenly as possible throughout the entire amount of grease forced out of the component being monitored.

The procedure for analyzing a grease sample with the Spectroil M depends upon whether the grease sample is relatively soft or if it is relatively stiff. If the grease is soft, simply fill the sample vessel with the grease. Perform the same analysis procedure as would be followed for a routine oil sample.

If the grease is too thick, the above method will not work because the grease will be so rigid that the rotating disc will pick up the sample vessel with the grease and will lift it off the sample table. In this case::

- 1) either heat the grease to a significantly higher temperature to soften the grease, or
- 2) smear the grease around the outside circumference of the disc, then fill a sample holder with base oil (0 ppm oil) and rotate the grease smeared disc through the cap filled with base oil following the routine analysis procedure

The analytical procedure for grease analysis consists of first analyzing a sample of a new, unused grease. Next analyze the used grease sample using either of the two methods described above. Abnormal wear will show up as an increase in the readings for Fe, Cu, Al, Pb, Sn, Si, etc., or whatever metals or additives are in the used grease.

3.0 OPERATOR MAINTENANCE

To maintain the Spectroil M performance, periodic maintenance must be performed by the operator. This maintenance falls into two categories:

- Daily Maintenance
- Scheduled Maintenance

3.1 Daily Operator Maintenance

This section details the maintenance actions required of the operator on a daily basis. These maintenance actions pertain mainly to the operator accessible assemblies such as the sample stand, the read-out and control panel, and the automatic printer. Any maintenance that is required to be performed in the excitation source, optics, or computer electronics is strictly limited to technically skilled personnel.

Table 3-1 lists each maintenance action that is authorized and required of the operator.

3.2 Scheduled Periodic Maintenance

This section of the manual details the maintenance action that is required to be performed on the Spectroil M at regularly scheduled intervals. Periodic maintenance will keep the instrument in good working condition and help to identify sources of future trouble before they cause serious downtime. The following tables separate periodic maintenance inspections by subassembly and item. Each item has instructions on what maintenance action is required, the interval in which the maintenance should be performed and by which maintenance level.

The scheduled maintenance tables are as follows:

- Table 3-2, External Housing Inspections
- Table 3-3, Internal Housing Maintenance Inspections
- Table 3-4, Excitation Source and Power Distribution Maintenance Inspections
- Table 3-5, Microprocessor Maintenance Inspections

TABLE 3-1
DAILY OPERATOR MAINTENANCE

CAUTION: Do not use alcohol or Chlorinated Solvents to clean plastic or painted surfaces.

Component	Required Maintenance	Frequency	Maintenance Level
Plate, Mounting, Sample Stand Component	Clean to remove oil and carbon buildup especially between disc electrode shaft and rod electrode holder. refer to Section 1.6.	Every 5 burns	Operator
Window, Quartz, Protective	Clean to remove oil and carbon splashes with isopropyl alcohol or an ammonia based window cleaner. Refer to Section 1.6.3.	Every 5 burns	Operator
Sensors, Sample Stand	Using a Q-tip, clean to remove oil and carbon splashes with isopropyl alcohol or an ammonia based window cleaner. Refer to Figure 1-4.	Daily	Operator
Sample Stand Area	Clean complete sample chamber to remove oil splashes and carbon buildup. Refer to Section 1.6.2.	Twice Daily	Operator
Door, Sample Stand	Clean complete door to remove oil splashes and carbon buildup. Refer to Section 1.6.	Twice Daily	Operator
Electrode Sharpener	Rotate cutting blade to new edge. (Can be performed until all three edges have been used.) Refer to Section 3.3.	As Required	Operator
Panel, Readout and Control	Inspect for oil splashes and carbon residue. If present, remove with mild cleaning detergent.	Daily	Operator
Frame and Exterior Panels	Inspect for oil splashes and dust buildup. If present, remove with mild detergent. CAUTION: DO NOT USE ALCOHOL OR CHLORINATED SOLVENTS TO CLEAN PLASTIC OR PAINTED SURFACES.	Daily	Operator
Printer	Inspect for worn ribbon, loose cable connectors, and dirt and dust buildup. Replace worn ribbon, tighten loose connections and clean accordingly. Refer to printer operation and maintenance manual.	Daily	Operator

TABLE 3-2
EXTERNAL HOUSING MAINTENANCE INSPECTIONS

Component	Required Maintenance	Frequency	Maintenance Level
Filter on Heat Exchanger	Inspect for dust and dirt buildup. Clean in detergent and water bath by swishing vigorously.	Weekly or as required depending on operating environment.	Operator
Filter, Sample Stand Exhaust	Inspect for dust and dirt buildup. Clean or replace if holes in the filter are blocked.	Weekly	Operator
Frame and Exterior Panels	Inspect for oil, dust, dents, scratches and rust. Clean with mild detergent and if necessary, sand and repaint.	Monthly	Operator
Hardware	Inspect for loose or missing hardware. Tighten loose hardware and replace rusted hardware.	Monthly	Operator
External Cables	Inspect for loose connections. Inspect for damage.	Monthly	Operator
Shaft, Disc Electrode	Clean residue (varnish) from splined end with an ink eraser.	Monthly	Operator

TABLE 3-3
INTERNAL HOUSING MAINTENANCE INSPECTIONS

Component	Required Maintenance	Frequency	Maintenance Level
Fans	Inspect for smooth rotation. Check for dust and dirt buildup on blades. Replace if binding is evident. Clean blades if necessary. Frequency - Six months	Six Months	Operator
Wiring	Inspect for broken or bent wiring connections. Inspect for frayed or burned insulation.	Six Months	Operator
Fuses	Inspect for open or over rated fuse usage. Replace as required.	Six Months or as Required	Operator
Transformers	Inspect for good electrical connection and signs of overheating.	Six Months or as Required	Operator
Video Monitor	Inspect for dust and dirt buildup. With a soft cloth, wipe clean if necessary.	Six Months	Operator
Signal Connectors	Inspect all connectors for proper seating in sockets. <i>CAUTION: Do not remove or connect any signal cables with the power on.</i>	Six Months	Operator

TABLE 3-4
EXCITATION SOURCE AND POWER DISTRIBUTION MAINTENANCE INSPECTION

Component	Required Maintenance	Frequency	Maintenance Level
Component Mounting Boards	Inspect the board for proper connector seating. Inspect for burn marks or discolored components. Inspect lower cabinet for dust and dirt buildup. Vacuum if necessary.	Six Months or as Required	Operator/Technician
Capacitors	Check each capacitor for signs of bulging, discolored containers or signs of leaking. Replace if signs of overheating are evident.	Six Months or as Required	Operator/Technician
Resistors	Check each resistor for signs of bulging or discoloration. Replace if overheating is evident.	Six Months or as Required	Operator/Technician
Contactors	Inspect for good electrical connection. Observe relay operation. Replace if intermittent or sluggish.	Six Months or as Required	Operator/Technician
Transformers	Inspect for signs of arc-over and overheating. Replace if evident.	Six Months	Operator/Technician
Analytical Gap	Inspect the rod electrode holder and gap setting device for smooth sliding and release. If tight or binding, adjust or remove and replace gap setting device, refer to the Maintenance Manual Section 2.4.1.1. Check the analytical gap distance. It should be 0.090 inches.	Six Months or Every 2,000 Burns	Operator/Technician
Auxiliary Gap	Polish tips to remove corrosion.	Six Months or Every 2,000 Burns	Operator/Technician
Auxiliary Gap	Check electrode shape. If electrode points are flat, remove electrodes and replace them. At sea level, reset the auxiliary gap distance to approximately 0.135 inches. Verify distance by checking excitation source frequency with Source Frequency Test Meter (refer to the Maintenance Manual Section 2.2.)	Six Months or Every 2,000 Burns	Operator/Technician
Motor, Disc Electrode	Check motor rotation. If loose or binding, align and tighten.	Six Months	Operator/Technician
Cables	Check high voltage cables in the excitation source for signs of arc-over or damage. Check electrical connections. Replace cables if arc-over is evident. Tighten connections if loose.	Six Months	Operator/Technician
Shaft, Disc Electrode	Check to be sure that the disc electrode shaft is tight and properly aligned. See the Maintenance Manual Section 2.4.1.2 for replacement procedure and tracking check.	Six Months	Operator/Technician

Disc Electrode to Rod Electrode Alignment	Check the alignment. If the rod point is more than 25% away from disc center, adjust. See the Maintenance Manual Section 2.4.1.3	Six Months	Operator/Technician
Auxiliary Gap Fan	Check for smooth rotation and cleanliness. Clean if necessary.	Six Months or Every 2,000 Burns	Operator/Technician

TABLE 3-5
MICROPROCESSOR MAINTENANCE INSPECTIONS

Component	Required Maintenance	Frequency	Maintenance Level
Printed Circuit Assemblies	Inspect each card for signs of discoloration due to component overheating. If present, replace the appropriate card. Check each card for proper connector seating. Re-seat if necessary.	Six Months	Operator/Technician
Optical Fibers	Inspect each fiber optic in the M58000 and M59200 cards for a tight mounting. If loose, remove the fiber and re-seat. If broken, replace. Frequency - Six months	Six Months	Operator/Technician
Cables, Interconnecting	Check both connectors of each interconnecting cable. If loose, Re-seat the connector in the appropriate plug. If broken or damaged, replace. CAUTION: Do not remove or connect any signal cables with the power on	Six Months	Operator/Technician

3.3 Procedure to Replace Electrode Sharpener Cutter Blade, M90102

The cutter blade has three sharpened edges and can be rotated and used three times before it has to be replaced.

To replace or rotate the cutter blade to a new cutting edge, unplug the sharpener power connector J2 at the power connection plate. Next, remove the graphite collector barrel assembly to empty out any accumulated graphite and to expose the cutter blade, Figure 3-1.

CAUTION: *The sharpener barrel will most likely be filled with graphite that has been removed from rod electrodes. Carefully follow the procedure below to avoid spilling graphite.*

To remove the barrel assembly, locate the sharpener over a waste basket with the collector barrel pointing downward. Grasp the barrel with the opposite hand and rotate it while pulling it away from the motor mount and face plate. Once the O-ring disengages the face plate, it will be easy to separate and empty.

Use a flat blade screw driver to remove the #4-40 screw which mounts the cutter blade, see Figure 2-14. When replacing or rotating the cutter blade, be sure to place the rear edge of the cutter blade tight against the cutter head. This is the reference point to achieve the correct angle on the graphite rod electrode. Replace the barrel assembly and reconnect the sharpener to its power connector.




Figure 3-1, Electrode Sharpener Blade

4.0 PERFORMING CALIBRATION CURVE VERIFICATION

The purpose of performing a calibration curve verification is to determine if the instrument repeats the curve generated at the factory or by an authorized service representative.

To perform a calibration curve verification, the instrument must first be standardized. Refer to Section 2.5 in this manual for a detailed procedure on daily standardization. When the instrument has been standardized using the calibration standards for the new program, the calibration curve verification can be performed.

The calibration curve verification consists of performing an analysis of each synthesized standard as if it were an unknown sample. It is recommended that the operator conduct ten analyses of each standard and perform statistics on the measurements to obtain the average and standard deviation for each element. Instrument performance for wear metal analysis should be within the limits listed in Tables 4-1 and 4-2.

4.1 Repeatability Testing

Perhaps one of the most important technical characteristics of a spectrometer is its ability to perform the same measurement over and over again with the same result. This characteristic is referred to as repeatability, reproducibility, sigma, standard deviation or precision. Repeatability is determined by the standard deviation of a series of measurements made on the same sample. Mathematically, standard deviation is calculated as:

$$\text{Std. Dev.} = \sqrt{[N(\sum X_i^2) - (\sum X_i)^2] / [N(N-1)]}$$

where:

N = the number of analyses (normally 10)

$\sum X_i^2$ = is the sum of the 10 squared individual measurements

$(\sum X_i)^2$ = is the square of the sum of the 10 individual measurements.

For the purpose of conveniently comparing the standard deviation to the mean for a series of measurements, relative standard deviation, or R.S.D., is used. R.S.D. expresses the standard deviation as a percent of the mean and is calculated as:

$$\text{R.S.D.} = 100[\text{Std. Dev.}/\bar{X}_{\text{avg}}]$$

The Spectroil M automatically calculates mean, standard deviation and R.S.D. for a series of measurements by pressing function key 5 (F5).

NOTE: At 0 ppm, R.S.D. is not considered a valid statistic. A minimum of three warm-up burns must be made prior to the performance of statistical analysis.

4.1.1 Repeatability Specifications

The Spectroil M is expected to perform within

repeatability specifications. The repeatability performance of the Spectroil M is part of the final test and calibration procedure done by Spectro before a Spectroil M/C or M/N is delivered. The specification to which the Spectroil M/C and M/N conform is summarized in Table 4-1. This table gives standard deviation values as a function of element and concentration. For example, if a repeatability test is conducted for titanium at 30 ppm, the standard deviation for ten measurements should be 1.87 ppm or less.

4.1.2 Repeatability Test

The level of repeatability given in Section 4.1.1 is obtained by burning the same standard ten times in succession. Press function key 5 (F5), left click the statistics icon, or select Operations/Statistics to obtain the mean, standard deviation and R.S.D. After five burns, the repeatability data can be checked by pressing the function key 5 (F5), left click the statistics icon, or select Operations/Statistics. If the repeatability is acceptable, burning five more times will, most probably, only

Table 4-1, Acceptable Repeatability Indices for Wear Metals - Standard Deviation

Conc. PPM	Fe, Al, Cr, Cu, Mg, Ni, Si	Ti, B	Ag, Na, Mo	Pb	Sn	Zn	V, Mn, Cd, Ca, Ba	P
0	0.50	0.50	0.50	0.90	1.00	0.50	0.5	N/A
5	0.56	0.58	0.64	0.95	1.04	0.78	N/A	N/A
10	0.71	0.78	0.94	1.08	1.17	1.30	0.75	N/A
30	1.58	1.87	2.45	2.01	2.06	3.63	1.75	5.5
50	2.55	3.04	4.03	3.13	3.16	6.02	2.75	7.5
100	5.03	6.02	8.02	6.07	6.08	12.0	3.75	10

Conc. PPM	Fe Ag Mo	Al Ni Si	Cr	Cu Mg	Na	Pb Sn Ti B	Zn	V	Mn	Cd Ba	Ca	P
300	24.0	15.0	27.0	48.0	18.0	36.0	10	12	10	10	15	
500	40.0	25.0	45.0	80.0	30.0	60.0	20	25	20	20	25	
700	56.0	35.0	63.0	112	42.0	84.0	N/A	N/A	N/A	N/A	N/A	
900	72.0	45.0	81.0	144	54.0	108	45	40	40	40	40	
5000	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	150	200
10000	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	300	300

improve the level of repeatability. (The nature of the standard deviation formula is such that the more tests that are done, the greater the divisor "N(N-1)" becomes, causing the standard deviation to become smaller.)

4.1.3 Factors Affecting Repeatability

In order to achieve this level of repeatability or better, the repeatability test must be done under ideal conditions. Many factors affect repeatability. Among them are:

1. The sample must be homogenous. The repeatability test done at Spectro is always done with standards. Routine samples are never used for repeatability testing because it cannot be assured they do not contain agglomerates, second phases and large particulates, all of which will affect repeatability.
2. The Spectroil M must be on profile. If analytical lines are off profile, the repeatability will be adversely affected. If the repeatability specifications cannot be met, one of the first diagnostic tests is to check profile (Section 2.7).
3. The quality and handling of the disc and rod electrodes will affect repeatability. The density, and hence the porosity, residual contamination and dimensional accuracy of the disc and rod electrodes will affect repeatability. Care must be taken to properly sharpen the rod electrode (Section 1.2). Proper care must also be exercised when installing the electrodes (Sections 1.3 and 1.4).
4. The sample must be homogenized by shaking before filling the sample holders.
5. The sample holders must be filled to the same level (Section 1.5).
6. Line voltage to the Spectroil M must be within specification.
7. Electronic stability of the Spectroil M will

affect repeatability.

8. Sample stand geometry will affect repeatability. The rod electrode to disc electrode gap distance, the quartz lens assembly to arc distance, the position of the fiber optic within the lens assembly mounting block, and the angle of the quartz lens assembly with respect to the arc will affect not only the intensity of the light entering the entrance slit of the polychromator (optical assembly), but will also affect the repeatability. The calibration of the Spectroil M at the factory optimizes these adjustments.
9. A variety of mechanical or electronic faults could degrade repeatability. Among these are faulty CCD chips, damage to the entrance slit, or damage to the fiber optic cable.

The operator has control over the first five factors. If care is taken to properly operate the Spectroil M and repeatability is still worse than specification and if the Spectroil M passes the Dark Current Test and is on profile, then it is recommended that Spectro Incorporated Field Service be consulted. It is strongly recommended that adjustments to the sample stand as described in 8 above be made only by Spectro personnel or at the direction of Spectro personnel.

4.2 Accuracy Testing

The Spectroil M is expected to perform within accuracy specifications in the same way that it performs within repeatability specifications. Accuracy is the ability of a spectrometer to give the correct concentration value of a standard. Table 4-2 gives acceptable accuracy readings for wear metal elements as a function of the concentration of the standard.

Column 1 of Table 4-2 gives concentration values in ppm. For example, if a 50 ppm multi-element standard is burned on the Spectroil M, the average of ten burns for Aluminum is expected to be

Table 4-2, Acceptable Accuracy Indices for Wear Metals - Mean

Conc, PPM	Al Cr Ni Si	Ti B	Fe Ag Mo	Cu Mg	Pb Sn	Zn	Na	V Mn Cd Ba	Ca	P
0	0.88	0.89	0.91	0.92	1.60	0.96	1.01	0.5	0.5	N/A
5	1.20	1.30	1.50	1.61	1.98	1.99	2.59	N/A	N/A	N/A
10	1.59	1.78	2.21	2.44	2.43	3.19	4.36	1.5	1.5	N/A
30	3.33	3.93	5.23	5.91	4.47	8.15	11.6	3.25	3.25	15.5
50	5.12	6.14	8.29	9.43	6.64	13.1	18.9	5.5	5.5	18.2
100	9.65	11.7	16.0	18.2	12.2	25.6	37.1	10	10	20
300	27.8	33.9	46.7	53.5	34.3	75.6	110	32	32	35
500	46.0	56.1	77.5	88.8	56.6	126	183	53	53	60
700	64.2	78.3	108	124	78.8	176	255	N/A	N/A	N/A
900	82.4	101	139	159	101	226	328	95	95	105
5000	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	500	500
10000	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	1000	1000

50 ppm, plus or minus 5.12 ppm. Therefore, an average in the range of 44.88 to 55.12 ppm would be acceptable. Zinc is a more difficult element to accurately measure. At 50 ppm, an acceptable average of ten burns is 50 ppm plus or minus 13.1 ppm. Therefore, an average in the range of 36.9 to 63.1 ppm would be acceptable.

The same factors that affect repeatability as described in Section 4.1.3 also affect accuracy. A complete standardization of the Spectroil M should be performed prior to testing for accuracy. The average of ten burns should be used to determine the reading of a particular standard. Accuracy failures at low concentrations may be due to contamination and spot impurities in the disc electrodes. The effect of the contamination and impurities can be compensated for by performing the electrode offset procedure in Section 3.3. Spectro Incorporated Field Service should be consulted if the Spectroil M is unable to meet the accuracy criteria presented in Table 4-2.

4.3 Disc Electrode Offset Procedure

The Spectroil M is designed to incorporate a background measurement and correction system. The purpose of this system is to offset or null the output of all CCD chips when measuring a 0 ppm standard. This is also known as measuring background light because 0 ppm has no concentration of elements present in the sample. Therefore, the light produced when analyzing a 0 ppm standard must only be background emission. This is, however, only true in theory.

In practice, elemental contamination is present in everything used for the analysis process. The sample holders may pick up contamination from the environment, the 0 ppm standard may have sub-ppm trace levels of certain elements, and the graphite disc electrodes are known to have trace contamination of certain elements. Manufacturers of graphite electrodes commonly list and quantify the known trace or spot impurities on each box of disc electrodes. The purpose of this procedure is to offset these trace contaminants in the consumables and is absolutely necessary to be performed.

This procedure should be performed every time a new batch and/or lot number of disc electrodes are to be used. For maximum efficiency in a laboratory operation, all graphite disc electrodes should be grouped and stored by batch and lot number. Only one batch or lot should be used at a time until it is totally consumed. Once a new lot is opened and the instrument is standardized to the new lot, the low end of the calibration curve (5 ppm and/or 10 ppm) should always be checked. If accuracy at these levels fails to meet the specified criteria, it may be due to variance in trace contaminants levels and the following procedure should be performed to correct for the presence of this contamination.

1. To perform the disc electrode offsets operation, the operator can select function key 10 (F10), left click the disc electrode offsets icon or choose Operations/Offsets/Perform Disc Offsets from the pull down menu options.

NOTE: The software has now been placed in a mode that is not to be used for normal operation. "Disc Electrode Offset" appears highlighted in red across the top of the analysis program screen to draw attention to the fact that this is not a normal condition for operation. Once this procedure is completed, the readout mode will have to be restored to the instrument's normal operating condition.

2. Prepare five sample holders of the 0 ppm oil standard. The screen will automatically clear any previous measurements when this mode is selected. Using the new batch or lot of disc electrodes, burn these five samples in accordance with the routine operating procedures outlined in Section 2.3.
3. After all five samples have been burned, press function key 6 (F6), the average icon, or Operations/Average to calculate the average of the five measurements. Once the average is calculated, a dialog with the title Background Correction Factors will appear

on the screen. This table will vary from instrument depending on the analytical configuration of the spectrometer; however, the format is the same as shown in Figure 4-1.

4. This dialog will display (from the left) ELEMENT in the second column, WAVELENGTH in the third, FORWARD intensity in the fourth, REVERSE intensity in the fifth, F/R RATIO in the sixth and BKG FACTOR in the last column. The cursor will appear in the upper right corner of the screen under the column BKG FACTOR. The absolute value for this mode is 1.00000 which indicates that the intensity produced in the forward or peak measurement is identical to the intensity produced in the reverse or background measurement.
5. The purpose of this procedure is to set the new ratio calculated and shown on the F/R RATIO column into the BKG FACTOR column for most elements.

NOTE: It is extremely important to pay attention to which elements this ratio is applied. Not all elements have a background factor as indicated by the value 0.00000. For the elements specified below, do not set the F/R RATIO value in the BKG FACTOR column. Failure to do so will adversely affect the analytical results of the instrument!*

Element	Wavelength	Forward	Reverse	F/R Ratio	Factor
1. Ag	328.068	0	0	0.00000	0.88049
2. Ag2	243.779	0	0	0.00000	0.00000
3. Al	308.216	0	0	0.00000	1.18456
4. B	249.678	0	0	0.00000	0.95434
5. Ba2	230.423	0	0	0.00000	0.00000
6. Ba	455.403	0	0	0.00000	0.91336
7. Ca2	448.478	0	0	0.00000	0.00000
8. Ca	393.365	0	0	0.00000	0.92830
9. Cd	226.502	0	0	0.00000	1.02476
10. Cr	425.435	0	0	0.00000	1.05217
11. Cu2	224.261	0	0	0.00000	1.11166
12. Fe	259.940	0	0	0.00000	0.98954
13. H	486.133	0	0	0.00000	0.00000
14. Mg	518.362	0	0	0.00000	1.10536
15. Mn	294.920	0	0	0.00000	1.04516
16. Mo	281.616	0	0	0.00000	1.12007

Figure 4-1, Sample Background Correction Factors Table

*Spectroil M - NaHi 568.861, P 510.656, Ca 445.478, and H 486.133

*Spectroil M/C & M/N - NaHi 568.861, P 510.656, Ca 445.478, MgHi 518.36, Ba 230.48, ZnHi 481.05 and H 486.133 and C 387.10.

*Spectroil M/F - MgHi 518.36, C 387.10 and H 486.133.

6. The cursor will automatically be located in the first row of the background factor column. To set the new F/R RATIO for all elements, left click the Set All button. Automatically, the new factors calculated and displayed in the F/R RATIO column will appear in the background factor column.
7. Once all of the elements that have a background factor have been updated to the new F/R RATIO, a copy of this screen may be made for your records by left clicking the Print button. Keep this printout for future reference. Left click the OK button to exit this dialog.
8. To exit the disc electrode offset procedure, press function key 10 (F10), left click the icon, or select Operations/Offsets/Perform Disc Offsets from the pull down menu options. This will place the software back to the Analysis Program screen which is the normal mode for operation and the mode the instrument was in before initiating this procedure.
9. This concludes the disc electrode offset procedure. Standardize the instrument for normal operation in accordance with Section 2.5 of this manual. Perform a daily standardization check in accordance with Section 2.4 to confirm that all elements at the lower concentrations meet the accuracy criteria. If any element fails to meet this criteria, contact Spectro Incorporated Field Service for assistance.

5.0 BACKUP OF COMPLETE OILM SOFTWARE

These steps will take the user through the process of backing up OilM for Windows on Microsoft® Windows XP.

1. While running OilMWindows open the File menu and select Exit. OilmWindows will shutdown.
2. Reopen OilMWindows by selecting from the Start menu or the shortcut on the desktop. After the communications are up select Operations/User Functions menu, Figure 4-2.
3. If user functions have been accessed, the window shown in Figure 4-5 open. Proceed with step 7. If no user functions have been previously accessed, the window as shown in Figure 4-3 will open. Proceed with step 4.
4. In the “Set up user functions paths” window click on the square button to the right of Removable drive:




Figure 4-2, Selection of User Function Menu




Figure 4-3, Set-up User Functions Paths Screen

5. The “Removable Drive:” window will open, Figure 4-4. Select the drive that says “Removable Disc (D:)”. With the disc drive highlighted click on the OK button.
6. The “Set up user functions paths” window, Figure 4-3, reopens. Click on the OK button to continue.
7. The “User Functions” window, Figure 4-5 opens. Select the “Backup to Removable Media” option. Click the OK button to continue.
8. OilMWindows will start a transfer of all OilM Windows and program files to the removable media.
9. When all the files are copied the “Backup completed successfully” message box will open. Click on the OK button to complete the back up procedure.




Figure 4-4, Removable Drive Selection Screen



Figure 4-5, Operations/User Functions Menu